

EDGE SPORT Assembly Instructions



3003 9th Ave. SW • Watertown, SD 57201
1-800-743-7738 • info@PivotalHealthSolutions.com

www.TeamEdgeAthletics.com



TOOLS NEEDED:

9/16" Wrench, 9/16" Socket and Phillips Head Screwdriver


IN THE BOX:



(1) Edge Sport Table




(1) H-Brace, (2) Leg Braces



(8) Bolts, (8) Nuts, (16) Washers, (4) Screws

Edge Sport Table Top



H-Brace

Leg Braces

TABLE PARTS:

EDGE SPORT Assembly Instructions



Unpack table from box. Turn table upside down on a protected surface.



Extend legs on both ends of the table into the upright position.



Install 2nd (lower) bolt on all (4) legs. Insert bolt, (2) washers and then nut. Hand tighten all bolts.



Install H-Brace into slit on the leg brace, with (4) 1-1/4" screws.



Tighten all (8) leg braces bolts using 9/16" socket and 9/16" wrench.



Turn table right side up.

Your Edge Sport table is now ready for use.

*NOTE: Routinely check table to make sure all hardware is tight and tighten accordingly when necessary.