

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



PHS Medical
by Pivotal Health Solutions

WWW.PHSMEDICALSOLUTIONS.COM

800-743-7738 • FAX: 605-882-8398

THERA-P™ ELECTRIC TREATMENT TABLE

STANDARD FEATURES:


- Load Capacity: 650 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9020/PT9028)
- White powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 27" Wide cushions in 7 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90° (PT1177: +35° to -45°)
- Lower extremity angle range: 0 to 65° (PT1177: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant

NOTE: International (240V) option available



Ease of use, multiple options and a sleek but durable design make the THERA-P line of electric treatment tables the choice of therapists everywhere. Seven top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The THERA-P electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system (with elevating foot pedal) and easy-to-operate caster system for better ergonomics.

TABLE CONFIGURATIONS:



☐ **PT1101**
1 Section
27" Width

Qty _____

☐ **PT1102**
2 Section
27" Width

Qty _____

☐ **PT1103**
3 Section
27" Width

Qty _____

☐ **PT1104**
4 Section
27" Width

Qty _____

☐ **PT1105**
5 Section
27" Width

Qty _____

☐ **PT1106**
6 Section
27" Width

Qty _____

☐ **PT1177**
7 Section
27" Width

Qty _____

THERA-P™ ELECTRIC TREATMENT TABLE

TABLE FEATURES:



OPTIONS:


- ☐ **PT9021** - 650 lb Lifting Capacity (240V for PT1101, PT1102, PT1103, PT1104, PT1105, & PT1106)
- ☐ **PT9029** - 650 lb Lifting Capacity (240V for PT1177)
- ☐ **PT303x** - 30"W Cushion Configuration
- ☐ **S1077** - Dual Foot Control Pedal
- ☐ **S1135** - Manual Elevating Midsection (Only available on PT1103 & PT1105)
- ☐ **EPTPH01** - Paper Roll Holder and Cutter (Only available on PT1101, PT1102, or PT1103)
- ☐ **S1078** - Dual Foot Bar Control
- ☐ **PT10001** - Nose Hole Plug



Dual Foot Control Pedal (x2) (S1077)



Manual Elevating Midsection (S1135)



Dual Foot Bar Control (S1078)



Paper Roll Holder and Cutter (EPTPH01)

STANDARD BASE COLOR:

- ☐  White

CUSTOM BASE COLORS:

- | | |
|---|---|
| ☐  Black | ☐  Green |
| ☐  Light Gray | ☐  Forest Green |
| ☐  Charcoal | ☐  Purple |
| ☐  Royal Blue | ☐  Yellow |
| ☐  Navy Blue | ☐  Orange |
| ☐  Red | ☐  Brown |
| ☐  Burgundy | ☐  Bright Blue |

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| ☐  Royal Blue | ☐  American Beauty |
| ☐  Imperial Blue | ☐  Lark |
| ☐  Blue Ridge | ☐  Burgundy |
| ☐  Deep Phoenician | ☐  Dove |
| ☐  Forest | ☐  Charcoal |
| ☐  Black | ☐  Mocha |

BRANDING KITS:

☐ **BRANDING KIT 1 - S1170**

- Slip Cover
- Custom Table Base Color
- Custom Table Base Graphics

**Not available on PT1101, PT1104, or PT1106*

☐ **BRANDING KIT 2 - S1171**

- Slip Cover
- White Table Base Color
- Custom Table Base Graphics

**Not available on PT1101, PT1104, or PT1106*

☐ **BRANDING KIT 3 - S1173**

- Deboss Logo
- Custom Table Base Color
- Custom Table Base Graphics

**Not available on PT1104 or PT1106*

☐ **BRANDING KIT 4 - S1174**

- Deboss Logo
- White Table Base Color
- Custom Table Base Graphics

**Not available on PT1104 or PT1106*

☐ **BRANDING KIT 5 - S1176**

- Color Edge Print
- Custom Table Base Color
- Custom Table Base Graphics

**Not available on PT1104 or PT1106*

☐ **BRANDING KIT 6 - S1177**

- Color Edge Print
- White Table Base Color
- Custom Table Base Graphics

**Not available on PT1104 or PT1106*

CUSTOM CUSHION BRANDING:

- ☐ **PT3027**
Logo Slip Cover
(27"W x 12"H)



- ☐ **PT3030**
Logo Slip Cover
(30"W x 12"H)

- ☐ **LOGO DB**
Deboss Logo
(12"W x 12"H)
DEBPLT12
Deboss Plate



- ☐ **LOGO CP**
Color Edge Print
(12"W x 12"H)



COLOR EDGE PRINT UPHOLSTERY CHOICES:

- | | |
|---|--|
| ☐  Blue Ribbon | ☐  Celestial |
| ☐  Purple Passion | ☐  Valentine |
| ☐  Crimson | ☐  Dolphin |
| ☐  Black | |

****Color Edge Print upholstery colors do not match
Standard Upholstery colors.**

CUSTOM BASE OPTIONS:

- ☐ **S1178**
Table Base Graphics
- ☐ **PT9020-CC**
Pick Custom Base Color
- ☐ **S2000**
Custom Table Base Powder Coat Color and Table Base Graphics



*PT1106 with Dual
Footbars and Table
Base Graphics*