

NAME: \_\_\_\_\_

PHONE: \_\_\_\_\_

SCHOOL/ORGANIZATION: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

DEALER/REP: \_\_\_\_\_

SALES REP: \_\_\_\_\_ QUOTE: \_\_\_\_\_



The Leg and Shoulder Treatment Cabinet supports bilateral treatment and specific anatomical positioning to achieve the desired results for patients, while the practitioner can store treatment equipment at arm's reach. The unique design of the CAB-100-LST allows for patient positioning and range of motion to perform a wide range of activities using one stable platform. To maximize the efforts of physical therapists, athletic trainers, and clinicians, this cabinet is a must for each rehabilitation practice or athletic training room.

**STANDARD FEATURES:**

- 30"W x 84"L x 34"H
- Load Capacity: 500 lbs.
- 3 Section cabinet
- (4) Hinged doors on each side
- Pass through open storage on foot end
- (1) Fixed height shelf in each section
- 70° Mechanical friction lock elevating back
- 90° Split leg elevation
- Drop out leg/shoulder area for improved range of motion
- Available in wood only
- 110° Euro hinges on doors
- 2" High density foam cushion
- Rounded corners
- Nose Hole
- Black wire pull handles
- US Patent: 10,857,052

*CAB-100-LST*  
**LEG & SHOULDER (LAST) TREATMENT CABINET**



**STANDARD CABINET:**

- CAB-100-LST**  
30"W x 84"L x 34"H  
Qty \_\_\_\_\_

**OPTIONS:**

- CAB-FL** - Flush Mount Handles on Doors  
 **Lock Double Door** - Lock on Double Doors



*Fixed height shelf*



*Manual elevating lift back*



**STANDARD WOOD CHOICES:**

-  Oak      Maple

**WOOD STAIN OPTION:**

- CAB-Stain

**STAIN CHOICES:**

-  Warm Copper      Rustic Fall  
  Deep Garnet      Warm Caramel  
  Fossil Gray      Black Beauty  
  Tungsten     Custom Wood Stain Color

**STANDARD UPHOLSTERY CHOICES:**

-  Royal Blue      American Beauty  
  Imperial Blue      Lark  
  Blue Ridge      Burgundy  
  Deep Violet      Dove  
  Forest      Charcoal  
  Black      Mocha  
  Mandarin Orange








**CUSTOM BRANDING:**

- Decal Custom Small**  
Small Logo Decal  
(6"W x 6"H) 
- Decal Custom Large**  
Large Logo Decal  
(12"W x 12"H) 

**CUSTOM CUSHION BRANDING:**

- LOGO COVER**  
Logo Slip Cover  
(Lift Back Only)  
(30"W x 12"H) 
- LOGO DB**  
Deboss Logo  
(12"W x 12"H)  
**DEBPLT12**  
Deboss Plate 
- LOGO CP**  
Color Edge Print  
(12"W x 12"H) 

**COLOR EDGE PRINT UPHOLSTERY CHOICES:**

-  Blue Ribbon      Celestial  
  Purple Passion      Valentine  
  Crimson      Dolphin  
  Black

*\*\*Color Edge Print upholstery colors do not match  
Standard Upholstery colors.*