

Name:	
School/Organization:	
Dealer/Rep:	


Phone:			
E-mail:			
Sales Rep:		Quote:	

PARALLEL BARS

PLATFORM MOUNT

STANDARD FEATURES:

- Load capacity: 500 lbs.
- Bars: 19-30"W and 26-44"H
- Platform 48"W
- 19-30"W x 84" or 120"L x 26-44"H
- Stainless steel hand rails (1.5" Diameter)
- Posi-Lok locking knobs



AM-710

STANDARD CONFIGURATIONS:

☐ AM-707 - 7' Platform Mount Parallel Bars 19-30"W x 84"L x 26-44"H - Platform 48"W	☐ AM-710 - 10' Platform Mount Parallel Bars 19-30"W x 120"L x 26-44"H - Platform 48"W 2 Piece Platform
QTY:	QTY:

OPTIONS:

☐ AM-714 - Abduction Board for Parallel Bars .75"W x 94"L x 5"H
QTY:
☐ AM-715 - Balance Beam for Parallel Bars 5"W x 80"L x 2.5"H
QTY:
☐ AM-716 - Foot Placement Ladder for Parallel Bars 24"W x 84"L x 4"H
QTY: