



SMART



A Pivotal Health Solutions Brand

stravasolutions.com

STRAVA SMART SPLINT BATH

A CLEAN, SMART, INNOVATIVE DESIGN

A state-of-the-art, first-of-its-kind splint bath precisely engineered to remove the guesswork and help prevent cross contamination

Automatic Dual Disinfection

The only splint bath designed with two built-in methods of disinfection to keep patients safe and meet regulatory compliance

Fast & Uniform Heating

Only 25 minutes to reach set water temperature with uniform and accurate heat distribution

Smart Technology

Preset and custom thermoplastic settings, auto on/off, low water sensor, alerts, water evacuation, and much more

kills
99.9%
of introduced
microbes

Meet Regulatory Compliance

PATENT PENDING





- ✓ UV-C Germicidal Kill
- ✓ High-Heat Sanitation Cycle
- ✓ Thermoplastic Presets
- ✓ Two-Stage Filtration
- ✓ Auto On/Off
- ✓ Fast & Uniform Heating
- ✓ Audible/Visual Alerts
- ✓ Bi-Fold Stainless Steel Lid
- ✓ Low Water Sensor
- ✓ Rapid Water Evacuation

THERMOPLASTIC PRESETS

Recommended time and temperature presets



Choose pre loaded thermoplastic settings from top manufacturers or create your own from the LCD touch screen – only available on the Strava Splint Bath.



DIGITAL TOUCH INTERFACE

Traditional splint bath controls are antiquated and imprecise

A full menu of features allow clinicians control and customization from an easy to use touch display



Fully retractable bi-fold lid fits under most cabinets for efficient use of counter top space

