

Name: \_\_\_\_\_

School/Organization: \_\_\_\_\_

Dealer/Rep: \_\_\_\_\_

Phone: \_\_\_\_\_

E-mail: \_\_\_\_\_

Sales Rep: \_\_\_\_\_ Quote: \_\_\_\_\_

# ALUMA ELITE

## CHANGING TABLE

### STANDARD FEATURES:

- Load capacity: 500 lbs.
- Welded aluminum frame
- Mounted hand control (only available on electric elevation)
- Drop down side rails
- H-brace construction
- 2" High density foam cushion
- Rounded corners
- Tough, anti-bacterial and anti-microbial powder coat finish
- Ships fully assembled
- Frame only available in Dusty Grey

### FIXED HEIGHT TABLE CONFIGURATIONS:

☐ **A4327-3052**  
31"W x 55"L x 27"H  
QTY: \_\_\_\_\_

☐ **A4327-3052ADA**  
31"W x 55"L x 22"H  
QTY: \_\_\_\_\_

### ELECTRIC HEIGHT TABLE CONFIGURATIONS:

☐ **A4327-3052EL**  
31"W x 55"L x 27-37"H  
QTY: \_\_\_\_\_

☐ **A4327-3052ELADA**  
31"W x 55"L x 22-32"H  
QTY: \_\_\_\_\_



A4327-3052EL

### OPTIONS:

- ☐ **AEBC-BS-B** - Bottom Shelf
- ☐ **A4725CT** - Sliding Door with Shelf
- ☐ **A4975** - Foot Switch
- ☐ **EPTPH01** - Paper Roll Holder (On Foot End With Paper Cutter on Head End)

### STANDARD UPHOLSTERY CHOICES:

- |                                       |                                       |
|---------------------------------------|---------------------------------------|
| <input type="radio"/> Royal Blue      | <input type="radio"/> American Beauty |
| <input type="radio"/> Imperial Blue   | <input type="radio"/> Lark            |
| <input type="radio"/> Blue Ridge      | <input type="radio"/> Burgundy        |
| <input type="radio"/> Deep Violet     | <input type="radio"/> Dove            |
| <input type="radio"/> Forest          | <input type="radio"/> Charcoal        |
| <input type="radio"/> Mandarin Orange | <input type="radio"/> Mocha           |
| <input type="radio"/> Black           | <input type="radio"/> Putty           |

