

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

AMBULATION TRAINING

RAMP AND CURB TRAINING SET STANDARD FEATURES:

- Load capacity: 350 lbs.
- Overall dimensions: 84"W x 44"W x 8"H
- Hardwood plywood construction
- Adjustable curbs store beneath platform
- Anti-slip safety tread on platform and all steps for increased safety

STANDARD SIZE:

- WCCT** - Ramp and Curb Training Set
Overall dimensions: 84"W x 44"W x 8"H
Curb 1: 32"L x 32"W x 2"H
Curb 2: 36"L x 32"W x 4"H
Curb 3: 40"L x 40"W x 6"H

QTY: _____



WCCT

ROLLING STEP STANDARD FEATURES:

- Load capacity: 300 lbs.
- Hardwood plywood construction
- 30"W x 30"L x 8"H
- Anti-slip safety tread on platform for increased safety
- (4) 2" In-line casters

STANDARD SIZE:

- WCRS** - Rolling Step
30"W x 30"L x 8"H

QTY: _____



WCRS

NESTED FOOT STOOLS STANDARD FEATURES:

- Load capacity: 350 lbs.
- Hardwood plywood construction
- Anti-slip safety tread on steps for increased safety

STANDARD SIZES:

- NFS-100** - Nested Foot Stool
Overall Dimensions: 16"W x 20"D x 8"H

Step 1: 10"W x 14"D x 2"H
Step 2: 12"W x 16"D x 4"H
Step 3: 14"W x 18"D x 6"H
Step 4: 16"W x 20"D x 8"H

QTY: _____

- NFS-101-16H** - Nested Foot Stool
Overall Dimensions: 16"W x 20"D x 16"H

QTY: _____



NFS-100



NFS-101-16H



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

STACKING STEP STOOLS STANDARD FEATURES:

- Load capacity: 500 lbs.
- Hardwood plywood construction
- Anti-slip safety tread on steps for increased safety

STANDARD SIZES:

- WCSS-100** - Stacking Step Stools
Step 1: 24"W x 12"D x 4"H
Step 2: 24"W x 12"D x 6"H
Step 3: 24"W x 12"D x 8"H
QTY: _____



WCSS-100

