

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

DOC

DECOMPRESSION TABLE

STANDARD FEATURES:

- Lifting capacity: 350 lbs.
- Electric elevation range: 21-29"
- Lateral flexion positioning and axial rotation
- Cervical capture headpiece
- Traction belting system enhances restraint for optimal decompression
- User-defined treatment hold and relax times
- Real-time digital treatment tracking
- Speed control
- Pre-programmed lumbar and cervical decompression protocols
- Digital Command Center includes:
 - Specific vertebral targeting
 - Separate lumbar and cervical traction/decompression programming
 - Continuous readout and graphing of treatment parameters
 - Interactive color display
- Reference library - including patient setup and training materials



E9011-B

STANDARD TABLES:

- E9011-B** - DOC Decompression Table (120V)
Table: 30"W x 77"L x 21"H
Space Requirement: 44"W x 102"L
QTY: _____
- E9014-B** - DOC Decompression Table (240V)
Table: 30"W x 77"L x 21"H
Space Requirement: 44"W x 102"L
QTY: _____
- E9022** - Tower Version (120V)
Table: 30"W x 77"L x 21"H
Tower: 27.5"W x 25"D x 72"H
Overall: 36"W x 144"L x 21"H
QTY: _____
- E9015** - Tower Version (240V)
Table: 30"W x 77"L x 21"H
Tower: 27.5"W x 25"D x 72"H
Overall: 36"W x 144"L x 21"H
QTY: _____

OPTIONS:

- E46500** - Elevation Foot Pedal
- E46510** - P-A Flexion Foot Pedal

STANDARD UPHOLSTERY CHOICES:

- | | |
|-------------------------------------|--------------------------------|
| <input type="radio"/> Royal Blue | <input type="radio"/> Lark |
| <input type="radio"/> Imperial Blue | <input type="radio"/> Burgundy |
| <input type="radio"/> Blue Ridge | <input type="radio"/> Dove |
| <input type="radio"/> Forest | <input type="radio"/> Charcoal |
| <input type="radio"/> Black | <input type="radio"/> Mocha |
| <input type="radio"/> Putty | |



E9022

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

DETAILS:

VIRTUAL COMMAND CENTER: 8" sleek, digital, color display slides to either side of the table.

DISC ANGLE LEVEL INDICATOR: Allows table to adjust to the exact angle for targeting a specific disc level increasing accuracy and specificity.

PATIENT SAFETY SWITCH: Safety feature allowing a patient to immediately stop the treatment due to pain, spasm or discomfort.


TABLE ELEVATION: A 21" to 29" range allows for customizable working height accommodating a variety of patient body types. The Americans with Disabilities Act allows this table to be eligible for the IRS tax code Section 44 credit in some cases. Check with your accountant for specific information.



CERVICAL CAPTURE HEADPIECE: Fully adjustable with the perfect combination of secure fit and patient comfort. Less discomfort, slipping, pressure than other decompression systems.



SECURE FIT PELVIC CAPTURE SYSTEM: The pelvic capture system securely captures the pelvis for a variety of body types providing just the right combination of fit and comfort.



Axial Rotation



Lateral Flexion



P-A Flexion



Knee & Cervical Bolster Set



Digital Command Center with Preprogrammed Options



Adjustable Arms



Adjustable Headpiece



"I looked at several spinal decompression tables before ultimately deciding on the DOC Decompression Table which was one of best business decisions I have made. The DOC Table is versatile, reliable and effective. Our patients love it. Results have been just as promised. The most exciting part is the table is making us money. It paid for it self in the first year. It is pretty much full all the time and we are just about to the point of needing a second table. I would recommend the DOC to anyone looking to add an effective, cash service to their practice."

Dr. Steve Brenan, Florida



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com