

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

EPD PACKAGE 2

EPIDURAL POSITIONING DEVICE & STAND

STANDARD FEATURES:

- Load capacity: 600 lbs.
- 20.5"W x 30"D x 50-60"H (Assembled)
- White powder coat base
- Black upholstery
- 180° Adjustable face rest
- 6 Position adjustable arm rest
- Arm rest size: 16.5"W x 10"D
- Arm rest adjustment: 32-43" from the floor
- Interchangeable torso cushions (1 flat and 1 elevated)
- Torso cushion size: 8"W x 8"H
- Torso cushion adjustment: 5.5-10"
- Positive locking clamp for EPD
- 7 Position adjustable foot rest
- Foot rest size: 20"W x 12"D
- Foot rest adjustment: 10.5-19.5" from the floor
- Medical grade profile locking casters
- Includes 1 Box (50) Sani-Cover® Fitted Disposable Face Rest Covers



EPD Package 2
(Epidural Positioning
Device + Stand)
with flat torso cushion



Elevated torso
cushion



FRC 30-50



Patents:
US 6,758,447,B2
and D455,495

STANDARD PACKAGE:

- EPD Package 2**
Epidural Positioning Device & EPD Stand
QTY: _____



EPD - Epidural
Positioning Device

LDS - EPD Stand

OPTIONS:

- FRC 30-50** Sani-Cover® Fitted Disposable Face Rest Covers (50 Covers per Box)
QTY: _____
- FRC 30-50-CASE** Sani-Cover® Fitted Disposable Face Rest Covers (50 Covers per Box) Case of 18 Boxes (900)
QTY: _____



Sani-Covers™ can be used
on the face rest and the arm
rest of the EPD

