

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

HT210 HAND THERAPY TABLE



STANDARD FEATURES:

- Load capacity: 264 lbs.
- Sit or stand
- Electrical elevation: 23-49”H
- 90 degree configuration
- Charcoal finish on base
- Anti-collision safety feature
- Hand Control with:
 - Constant smooth movement and soft start/stop
 - (3) Programmable memory locations
 - (2) USB ports
- Pre-assembled frame kits significantly reduces installation time
- (6) Leveling feet
- Some assembly required



HT210-7229L

TABLE CONFIGURATIONS:

- HT210-7229L**
72”L x 29”D x 23-49”H
Overall Footprint: 72”W x 72”D
QTY: _____


- HT210-7229R**
72”L x 29”D x 23-49”H
Overall Footprint: 72”W x 72”D
QTY: _____

OPTIONS:

- RE410** - Battery

TOP COLORS:

-  Brushed Concrete
-  Weathered Grey
-  Mineral Grey
-  Briarsmoke
-  Cocoa Bean
-  Nightshade
-  Sandlewood
-  White Wood Grain



Add one of each for a fully encompassing set up perfect for larger rooms



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com