

Name: _____

School/Organization: _____

Dealer/Rep: _____

Phone: _____

E-mail: _____

Sales Rep: _____ Quote: _____

PRO 6 LEG CLASSIC

WOOD TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 800 lbs.
- Solid hardwood 6 leg h-brace construction
- 2" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- Bottom shelf for added storage
- 70° Manual lift back (lift back option only)
- 6 Position elevation
- 45° Split leg elevation



FESS-6L-PRO with
Custom Stain

FLAT TOP TABLES:

- ☐ **FESS-6L-PRO**
Flat Cushion
30"W x 84"L x 31"H
QTY: _____

- ☐ **FESS-6LS-PRO**
Flat, Split Leg Cushion
30"W x 84"L x 31"H
QTY: _____

LIFT BACK TABLES:

- ☐ **TESS-6L-PRO**
Lift Back Cushion
30"W x 84"L x 31"H
QTY: _____

- ☐ **TESS-6LS-PRO**
Lift Back, Split Leg Cushion
30"W x 84"L x 31"H
QTY: _____

Manual Lift Back with
6 Position Elevation



OPTIONS:

- ☐ **ESS Cabinet 6L** - Shelf and Open Cabinet
☐ Head End ☐ Foot End
- ☐ **ESS Drawers 6L** - Shelf, Cabinet & (2) Drawers with Black Wire Pull Handles
☐ Head End ☐ Foot End
- ☐ **EWT-NH** - Nose Hole in Cushion
- ☐ **TSFO-FHDR** - Flush Mount Handles on Drawers
- ☐ **EPTPH01** - Paper Roll Holder (On Foot End with Paper Cutter on Head End)



Shelf, Cabinet & (2)
Drawers with Black
Wire Pull Handles

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

OPTIONS:

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|---|
| ☐  Royal Blue | ☐  American Beauty |
| ☐  Imperial Blue | ☐  Lark |
| ☐  Blue Ridge | ☐  Burgundy |
| ☐  Deep Violet | ☐  Dove |
| ☐  Forest | ☐  Charcoal |
| ☐  Mandarin Orange | ☐  Mocha |
| ☐  Black | ☐  Putty |

STANDARD WOOD CHOICES:

- | |
|---|
| ☐  Oak |
| ☐  Maple |

WOOD STAIN OPTIONS:

- ☐ **ESS-Stain**

CUSTOM STAIN CHOICES:

- | | |
|---|--|
| ☐  Warm Copper | ☐  Rustic Fall |
| ☐  Deep Garnet | ☐  Warm Caramel |
| ☐  Fossil Gray | ☐  Black Beauty |
| ☐  Tungsten | ☐ Custom Wood Stain Color: _____ |

CUSTOM BRANDING:

- ☐ **Decal Custom Small**
Small Logo Decal
(6"W x 6"H)



- ☐ **Decal Custom Large**
Large Logo Decal
(12"W x 12"H)



CUSTOM CUSHION BRANDING:

- ☐ **TIS 30 COVER/LOGO**
Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)



- ☐ **Logo DB**
Deboss Logo
(12"W x 12"H)
- ☐ **DEBPLT12 - Deboss Plate**



- ☐ **Logo CP - Color Edge Print**
(12"W x 12"H)



COLOR EDGE UPHOLSTERY CHOICES:

- | | |
|--|---|
| ☐  Blue Ribbon | ☐  Celestial |
| ☐  Purple Passion | ☐  Valentine |
| ☐  Crimson | ☐  Dolphin |
| ☐  Black | |

****Color Edge Print upholstery colors do not match Standard Upholstery colors.**