

Name: \_\_\_\_\_

Phone: \_\_\_\_\_

School/Organization: \_\_\_\_\_

E-mail: \_\_\_\_\_

Dealer/Rep: \_\_\_\_\_

Sales Rep: \_\_\_\_\_ Quote: \_\_\_\_\_

# THERA-P™ BARIATRIC ELECTRIC TREATMENT TABLE

## STANDARD FEATURES:


- Load capacity: 650 lbs.
- 34" or 36"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9024/PT9025)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 34" or 36" Wide cushions in 4 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant

**NOTE:** International (240V) option available



PT3613 with Dual Foot Bars

## TABLE CONFIGURATIONS:



1 Section    2 Section    3 Section    5 Section

**PT3011**  
1 Section  
34" Width

QTY: \_\_\_\_\_

**PT3012**  
2 Section  
34" Width

QTY: \_\_\_\_\_

**PT3013**  
3 Section  
34" Width

QTY: \_\_\_\_\_

**PT3015**  
5 Section  
34" Width

QTY: \_\_\_\_\_

**PT3611**  
1 Section  
36" Width

QTY: \_\_\_\_\_

**PT3612**  
2 Section  
36" Width

QTY: \_\_\_\_\_

**PT3613**  
3 Section  
36" Width

QTY: \_\_\_\_\_

**PT3615**  
5 Section  
36" Width

QTY: \_\_\_\_\_



Name: \_\_\_\_\_

Phone: \_\_\_\_\_

School/Organization: \_\_\_\_\_

E-mail: \_\_\_\_\_

Dealer/Rep: \_\_\_\_\_

Sales Rep: \_\_\_\_\_ Quote: \_\_\_\_\_

# OPTIONS:

## TABLE FEATURES:

- PT9025** - 650 lb Lifting Capacity
- S1136** - Manual Elevating Midsection (only available on PT3013, PT3015, PT3613, & PT3615)
- S1077** - Dual Foot Control Pedal
- S1078** - Dual Foot Bar Pedal



Manual Elevating Midsection (S1136)



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)

## STANDARD FRAME COLORS:

-  White

## CUSTOM BASE COLORS:

- |   |  |
|---|--|
| <input type="radio"/>  Black       | <input type="radio"/>  Green        |
| <input type="radio"/>  Light Gray  | <input type="radio"/>  Forest Green |
| <input type="radio"/>  Charcoal    | <input type="radio"/>  Purple       |
| <input type="radio"/>  Bright Blue | <input type="radio"/>  Yellow       |
| <input type="radio"/>  Royal Blue  | <input type="radio"/>  Orange       |
| <input type="radio"/>  Navy Blue   | <input type="radio"/>  Brown        |
| <input type="radio"/>  Red         | <input type="radio"/>  Sky Blue     |
| <input type="radio"/>  Burgundy    |  |

## STANDARD UPHOLSTERY CHOICES:

- |   |  |
|---|--|
| <input type="radio"/>  Royal Blue        | <input type="radio"/>  American Beauty |
| <input type="radio"/>  Imperial Blue     | <input type="radio"/>  Lark            |
| <input type="radio"/>  Blue Ridge       | <input type="radio"/>  Burgundy       |
| <input type="radio"/>  Deep Violet     | <input type="radio"/>  Dove          |
| <input type="radio"/>  Forest          | <input type="radio"/>  Charcoal      |
| <input type="radio"/>  Mandarin Orange | <input type="radio"/>  Mocha         |
| <input type="radio"/>  Black           | <input type="radio"/>  Putty         |

## CUSTOM BASE OPTIONS:

- S1178B** - Table Base Graphics
- PT9024-CC** - Pick Custom Base Color
- S2000B** - Custom Table Base Powder Coat Color and Table Base Graphics

