

Name: _____
 School/Organization: _____
 Dealer/Rep: _____

Phone: _____
 E-mail: _____
 Sales Rep: _____ Quote: _____

THERA-P™

ELECTRIC TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 650 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9020/PT9028)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 27"W Cushions in 8 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90° (PT1177: +35° to -45°)
- Lower extremity angle range: 0 to 65° (PT1177: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant

NOTE: International (240V) option available

TABLE CONFIGURATIONS WITH MANUAL LIFT BACK:





						
☐ PT1101 1 Section 27" Width QTY: _____	☐ PT1102 2 Section 27" Width QTY: _____	☐ PT1103 3 Section 27" Width QTY: _____	☐ PT1104 4 Section 27" Width QTY: _____	☐ PT1105 5 Section 27" Width QTY: _____	☐ PT1106 6 Section 27" Width QTY: _____	☐ PT1177 7 Section 27" Width QTY: _____
			☐ PT1104-SL 4 Section with Split Leg 30" Width QTY: _____			

TABLE CONFIGURATIONS WITH ELECTRIC LIFT BACK:

		
☐ PT1102-EL 2 Section 27" Width QTY: _____	☐ PT1103-EL 3 Section 27" Width QTY: _____	☐ PT1105-EL 5 Section 27" Width QTY: _____



PT1104-SL with Dual Foot Bar
Control and Branding Kit 3



STANDARD UPHOLSTERY CHOICES:

☐ Royal Blue	☐ American Beauty
☐ Imperial Blue	☐ Lark
☐ Blue Ridge	☐ Burgundy
☐ Deep Violet	☐ Dove
☐ Forest	☐ Charcoal
☐ Mandarin Orange	☐ Mocha
☐ Black	☐ Putty




OPTIONS:

TABLE FEATURES:

 Flat Cushion	 Lower Extremity Angle Ranges from 0 to 65°	 Lower Extremity Angle Ranges from 0 to 65° Optional Elevating Midsection Head Section Angle Ranges from +40° to -90°	 Lower Extremity Angle Ranges from 0 to 65° Head Section Angle Ranges from +40° to -90°
PT1101	PT1102/PT1102-EL	PT1103/PT1103-EL	PT1104
 *Table shown with options	 Lower Extremity Angle Ranges from 0 to 65° Each Armrest Vertical Ranges from 0 to 7° Head Section Angle Ranges from +40° to -90° Optional Elevating Midsection Each Armrest Horizontal Ranges from 0 to 45°	 Lower Extremity Angle Ranges from 0 to 65° Each Armrest Vertical Ranges from 0 to 7° Head Section Angle Ranges from +40° to -90° Each Armrest Horizontal Ranges from 0 to 45°	 Lower Extremity Angle Ranges from 0 to 75° Standard Electric Elevating Midsection Head Section Angle Ranges from +35° to -45° Adjustable Armrests Shoulder Dropouts 0 to -90°
PT1104-SL	PT1105/PT1105-EL	PT1106	PT1177

OPTIONS:

☐ PT9021 - 650 lb Lifting Capacity (240V FOR PT1101, PT1102, PT1103, PT1104, PT1105, PT1106)
☐ PT9029 - 650 lb Lifting Capacity (240V FOR PT1177)
☐ PT303X - 30"W Cushion Configuration (not available on PT1177)
☐ S1135 - Manual Elevating Midsection (only available on PT1103 & PT1105)
☐ S1077 - Dual Foot Control Pedal
☐ S1078 - Dual Foot Bar Control
☐ S2078 - Dual Foot Bar Control With Front Foot Bar Adder
☐ S1151 - (1) Battery & (1) Wall Charger
☐ S1152 - (2) Batteries & (1) Wall Charger
☐ EPTPH01 - Paper Roll Holder (only available on PT1101, PT1102, or PT1103)
☐ PT1001 - Nose Hole Plug

 Manual Elevating Midsection (S1135)	 Dual Foot Control Pedal (x2) (S1077)	 Dual Foot Bar Control (S1078)
 Dual Foot Bar Pedal With Front Foot Bar Adder (S2078)	 (1) Battery & (1) Wall Charger (S1151)	 Paper Roll Holder (EPTPH01)

BRANDING KITS:

☐ BRANDING KIT 1 - S1170 - Slip Cover - Custom Table Base Color - Custom Table Base Graphics	☐ BRANDING KIT 4 - S1174 - Deboss Logo - White Table Base Color - Custom Table Base Graphics
☐ BRANDING KIT 2 - S1171 - Slip Cover - White Table Base Color - Custom Table Base Graphics	☐ BRANDING KIT 5 - S1175 - Color Edge Print - Custom Table Base Color - Custom Table Base Graphics
☐ BRANDING KIT 3 - S1173 - Deboss Logo - Custom Table Base Color - Custom Table Base Graphics	☐ BRANDING KIT 6 - S1176 - Color Edge Print - White Table Base Color - Custom Table Base Graphics

CUSTOM BASE OPTIONS:

☐ S1178 - Table Base Graphics
☐ PT9020-CC - Pick Custom Base Color
☐ S2000 - Custom Table Base Powder Coat Color and Table Base Graphics



STANDARD BASE COLORS:

☐  White

CUSTOM BASE COLORS:

☐  Black	☐  Burgundy
☐  Light Gray	☐  Green
☐  Charcoal	☐  Forest Green
☐  Bright Blue	☐  Purple
☐  Royal Blue	☐  Yellow
☐  Navy Blue	☐  Orange
☐  Red	☐  Brown
☐  Sky Blue	

CUSTOM CUSHION BRANDING:

☐ PT3027 - Logo Slip Cover (27"W x 12"H)		☐ LOGO DB - Deboss Logo (12"W x 12"H)	
☐ PT3030 - Logo Slip Cover (30"W x 12"H)		DEBPLT12 - Deboss Plate	
☐ A4451 - Logo Slip Cover (PT1104-SL) (30"W x 12"H)			
☐ LOGO CP - Color Edge Print (12"W x 12"H)			

COLOR EDGE UPHOLSTERY CHOICES:

☐  Blue Ribbon	☐  Celestial
☐  Purple Passion	☐  Valentine
☐  Crimson	☐  Dolphin
☐  Black	

**Color Edge Print upholstery colors do not match Standard Upholstery colors.