

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

THERA-P™

TRACTION TABLE

STANDARD FEATURES:

- Load capacity: 650 lbs.
- 27"W x 78"L x 20-37"H
- Durable and stylish steel white powder coat base
- Whisper quiet lift system with dual foot bar control
- 27" Wide cushions in 2 different top configuration options
- 2" High density foam cushions
- Rounded corners
- Friction-free gas shock cushion separation with lock and unlock level
- (2) Low profile locking casters
- Head section angle range: +30° to -90°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- Lift back: 0 to 70°
- Table comes standard with traction head arm and mounting plate



TXT-1004

STANDARD TABLES:



TXT-1004

4 Section

Table Dimension: 27"W x 78"L x 20-37"H

Overall Dimension: 110"W x 64"W x 20-37"H

QTY: _____



TXT-1006

6 Section

Table Dimension: 27"W x 78"L x 20-37"H

Overall Dimension: 110"W x 64"W x 20-37"H

QTY: _____



4 Section



6 Section

Name: _____

School/Organization: _____

Dealer/Rep: _____

Phone: _____

E-mail: _____

Sales Rep: _____ Quote: _____

OPTIONS:

OPTIONS:

- ☐ **TXT-ACC** - Accessories Kit
Traction Head
Cervical Attachment
Traction Belt
Knee and Cervical Bolster Set
(1) Cervical Head Pillow
- ☐ **S1354** - Strapping System
Traction Belt
Knee and Cervical Bolster Set
(1) Cervical Head Pillow
- ☐ **S1353** - Traction Belt Only

TRACTION HEAD FEATURES:

- Load capacity: 650 lbs.
- Crystal clear displays
- Static, Intermittent, and Cyclic Traction
- Variable speed control
- Minimum and maximum force settings 2.5kgs/5.5lbs - 91kgs/200lbs
- Patient safety switch with audible alarm
- Traction force compensation
- Audible alarm protection monitoring of traction forces
- Continuous electronic monitoring of traction forces



Knee & Cervical
Bolster Set and
Traction Belt



Traction Head



Cervical
Attachment &
Head Pillow

STANDARD FRAME COLORS:

- ☐ White

CUSTOM BASE COLORS:

- | | |
|-----------------------------------|------------------------------------|
| ☐ Black | ☐ Green |
| ☐ Light Gray | ☐ Forest Green |
| ☐ Charcoal | ☐ Purple |
| ☐ Bright Blue | ☐ Yellow |
| ☐ Royal Blue | ☐ Orange |
| ☐ Navy Blue | ☐ Brown |
| ☐ Red | ☐ Sky Blue |
| ☐ Burgundy | |

STANDARD UPHOLSTERY CHOICES:

- | | |
|---------------------------------------|---------------------------------------|
| ☐ Royal Blue | ☐ American Beauty |
| ☐ Imperial Blue | ☐ Lark |
| ☐ Blue Ridge | ☐ Burgundy |
| ☐ Deep Violet | ☐ Dove |
| ☐ Forest | ☐ Charcoal |
| ☐ Mandarin Orange | ☐ Mocha |
| ☐ Black | ☐ Putty |

BASE OPTIONS:

- ☐ **S1178** - Table Base Graphics
- ☐ **PT9020-CC** - Pick Custom Base Color
- ☐ **S2000** - Custom Table Base Powder Coat Color and Table Base Graphics

