Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

PORTABLE
PARALLEL BARS

STANDARD FEATURES:

- Load capacity: 600 lbs on the base
  - 300 lbs per bar
- Foldable with casters for easy transport
- Transport mode:
  - 12.25-17.5"W x 60"L x 39.5"H
  - 360° Turn radius
- Extended mode:
  - 45.5"W x 60"L x 28.5-41"H
  - Adjustable handrail: 22.75-26.5"W
- Dual lock mechanism for height adjustment
- (4) 4"Total locking medical grade caster system
- Aluminum handrails
- Welded aluminum frame and walking surface
- Durable powder coat
- Anti-slip platform
- No assembly required



PTPB1000

STANDARD CONFIGURATIONS:

◯ PTPB1000 - Portable Parallel Bars 12.25 - 17.5"W x 60"L x 39.5"H
QTY:

Watch Set Up and Take Down Video, Get User Manual, and More...



Armedica
BY PIVOTAL HEALTH SOLUTIONS

800.701.5122
armedicamfg.com

Rev 02.20.2025